



## SPECIFICATIONS

Automac filling and sealing machine for glass ampoules "ABF-01"

Modifications: 1-2ml, 5-10ml, 20ml

Productivity (1-2 ml): 2000 - 2400 ampoules per hour.

Productivity (5-10 ml): 1000 - 1200 ampoules per hour.

Productivity (20 ml): 600 - 800 ampoules per hour.

Power: 0.37 kW

Gas consumption: 2-2.5 m<sup>3</sup> per hour

Air consumption: 0,08-0,1 MPa

Supporting combustion gas consumption: 0,7 m<sup>3</sup> per hour

Oxygen consumption: 0.1 MPa

Power supply: 220V or 110V,

Overall dimensions: 1300mm x 700mm x 1400mm

Weight: 150 kg

Weight with packing: 220 kg

 [www.minipress.ru](http://www.minipress.ru)

## ADDITIONAL INFORMATION

A unique series of automac machines for filling and sealing glass ampoules. Models differ from each other with small options, dosage and performance. The sampling of the ampoules comes from the hopper, then they get on the transporter in the filling. The operator adjusts the dosage, within the capabilities of the machine. Filling in glass ampoules is carried out by the group. The machine brings a group of ampoules under the needles, the conveyor stops and the needles fall into ampoules. Next comes the filling of ampoules. Then the glass ampoules are moved to the burner zone. The combustion temperature and the intensity of the flame are also regulated. After the neck of the ampoule has been heated by special clamps, the apex of the ampoule is formed. The tip has a rounded contour of regular uniform shape. The height of all ready-made ampoules and roundness is repeated constantly, without discontinuity. Therefore, the machine for filling and sealing glass ampoules requires a machine for washing and drying ampoules. And also for the internal injection machine for sterilization.

**PRODUCTION TIME**

**30-40 DAYS**

**ESTIMATED DELIVERY TIME**

**30-40 DAYS**

**PRICE**

**6779,66 USD**

**VAT**

**1220,34 USD**

**VALUE 8000.00 USD**

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank

**OUR SERVICES FOR CLIENTS:**

**1) TRAINING TECHNOLOGIES OF BOTTLING AND BATCHING AMPUL.**

We training on our equipment before buying. Tesng of materials. We provide consultaons on lling and sealing of ampoules, with a full range of addional equipment and consumables.

**2) WE COMPLETE PHARMACEUTICAL AND CONFECTIONERY PRODUCTION.**

We have automac and desktop machines for lling and sealing ampoules In Russia, which we provide to our customers. You can always contact us with your ideas and requests.

**3) SUPPLY OF EQUIPMENT FOR BOTTLING AND STORING OF GLASS AMPOULS.**

The products fully comply with the new GMP requirements. The release of equipment for sealing glass ampoules is accompanied by strict quality control, vericaon of all documentaon, well-organized aer-sales service of equipment.

**4) SPARE PARTS SUPPLY**

The manufacturer can also oer models of automac machines for lling and sealing ampoules, based on the goals and objecves of the customer. We will assist in the maintenance and repair of any equipment for sealing ampoules. Service center in Russia for the CIS countries.



*Roman Tsibulsky*

Owner and supervisor

**ROMAN TSIBULSKY**