

Commercial proposal

Compiled:	15.02.2018
Valid:	12 month
Author:	Roman Tsibulsky

Granulator for wet granulation of powders ZL-10

DESCRIPTION

Laboratory granulator of desktop type for wet granulation of powder masses in pharmaceutical manufacture. Manufacture up to 10 kg per hour. The sizes of ready granules are 0.8-1.5 mm. Material - stainless steel. Ideal for research laboratories and pilot production. We provide a full range of services: installation, training, start-up, repair. Simple mechanical design. Instruction for use in English and Russian.

Complies with the GMP standard. We provide detailed instructions for setting up this granulator model. Prior to shipment, granulators are tested and tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. The price is indicated taking into account customs payments in Russia and delivery to the city of the client.

MANUFACTURER - CHINA



MINIPRESS.RU



MINIPRESS.RU

SPECIFICATIONS

Granulator for wet granulation of powders "ZL-10"

Performance

Wet granulation: (with a screen size of 1 mm) 1 to 10 kg per hour

Replacement sieves: stainless

Dimensions of granules: from 0.5 mm to 3 mm

Accommodation: table

Engine power: 0,3 kW, 220V, 50 Hz

Rotational speed of the blade: 46 rpm

Overall dimensions: 370 mm x 320 mm x 570 mm

Weight: 20 kg

Shipping Weight: 30 kg

 www.minipress.ru

ADDITIONAL INFORMATION

Tabletop design of the granulator. At the base of the machine, the shaft of the electric motor on which the blades are clad. The walls of the casing into which the moistened powder is fed are made of perforated steel. As the blades rotate, the wet powder is wiped through the holes, the resulting granules can have both a sphere shape and worms. The shape of the granules depends on the properties of the powder, humectant and binder. Be sure to consult with us on the technology of preparing the powder for granulation. Granulator: does not mix powder; does not moisturize the powder; Do not mix the moistened powder; Do not dry the granules; Do not separate granules into good and substandard; Do not produce pellets of the same size shape. The function of the granulator is to wipe the moistened powder through the holes. Check whether the granulator is suitable for you is easy.

PRODUCTION TIME

30-40 DAYS

PRICE

13559,32 USD

ESTIMATED DELIVERY TIME

30-40 DAYS

VAT

2440,68 USD

OUR SERVICES FOR CLIENTS:

1) TRAINING GRANULATION TECHNOLOGIES.

You can train on our equipment before buying. Testing of materials. We provide consultations on the organization of granulation of materials, with a full range of additional equipment and consumables.

2) WE COMPLETE LABORATORIES FOR PHARMACEUTICAL PRODUCTION.

We have a laboratory test granulator for the production of pilot batches in Russia, which we provide to our customers. You can always contact us with your ideas and requests. Currently, the manufacturer of granulators can offer equipment for large-scale production of granules by all means.

3) DELIVERY OF GRANULATORS.

The products fully comply with the new GMP requirements. The release of granulators is accompanied by strict quality control, verification of all documentation, well-organized after-sales service of equipment.

4) SPARE PARTS SUPPLY

The manufacturer can also offer a granulator model based on the customer's goals and objectives. We will assist in the maintenance and repair of any equipment for granulation. Service center in Russia for the CIS countries.

VALUE 16000.00 USD

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank



Roman Tsibulsky

Owner and supervisor

ROMAN TSIBULSKY