

Commercial proposal

Compiled:	15.02.2022
Valid:	12 month
Author:	Roman Tsibulsky

Analyzer of dissolution of tablets and capsules RC-01

DESCRIPTION

Dissolution Analyzer is used to measure the rate and degree of dissolution of tablets, capsules, etc. under laboratory conditions. It has one tank and one blade located in one row. The main parts of this equipment rotate smoothly and reliably. The reservoir and blades are made of stainless steel (SUS316L). The water flow is circulated by a magnetic pump, distributes it evenly throughout the system, and the "water bath" maintains the same temperature.

This device for dissolving tablets is equipped with automatic temperature control, changes in speed and more. We provide detailed instructions for adjusting the dissolution analyzer. Before dispatch, dissolution analyzers are tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. The price is indicated taking into account customs payments in Russia and delivery to the city of the client.

MANUFACTURER - CHINA



MINIPRESS.RU



MINIPRESS.RU

SPECIFICATIONS

Dissolution analyzer for tablets and capsules "RC-01"
Speed range: (20-200) rpm
Accuracy of rotation speed: ± 2 rpm
Temperature of thermostat liquid: approximately $+ 45$ °C
Temperature stability: ± 0.3 °C
Time at which the mer is set: (5-900) min
Blade stopping time: ± 0.5 min
Time to stop the tank: ± 1.0 min
Power: 220V / 50Hz / 600W or 110V / 60Hz / 600W
Dimensions: 260mm x 300mm x 480mm
Size with packing: 520mm x 490mm x 565mm
Weight: 15 kg
Weight with packing: 20 kg

 www.minipress.ru

ADDITIONAL INFORMATION

The main parts of this equipment rotate smoothly and reliably. The reservoir and blades are made of stainless steel (SUS316L). The water flow is circulated by a magnetic pump, distributes it evenly throughout the system, and the "water bath" maintains the same temperature. This device for dissolving tablets is equipped with automatic temperature control, changes in speed and time. These parameters are set by the user. They can be controlled using the buttons. The results are displayed in the panel. The set values on the user and the indicators during the operation of the device can be displayed alternately. Automats: auto-test, auto-diagnosis, auto-signal. The equipment can test itself and the signal automatically. The dissolution parameters of the substances are stored in the equipment memory and can be displayed at any time. Time can be set on the user for many positions. The alarm device automatically sends a signal to any position of the user and then the tester automatically moves to the next position.

PRODUCTION TIME	30-40 DAYS
ESTIMATED DELIVERY TIME	30-40 DAYS

PRICE	2372,88 USD
VAT	427,12 USD

VALUE 2800.00 USD

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank

OUR SERVICES FOR CLIENTS:

1) TRAINING OF MEDICINE VERIFICATION TECHNOLOGY.

To all our customers who purchased equipment, we provide literature on the basics of encapsulation of powders. We carry out consultations by phone and by correspondence, we help to solve arising difficulties at manufacture on our equipment. We supply spare parts and consumables. We render services in the help on realization of the capsular equipment used.

2) WE COMPLETE PRODUCTION CAPSULES.

We have a wide range of additional equipment for capsule production. Upon customer's request, we will select a model of capsule filling equipment, granulators, powder mixers, dust removal equipment and capsule polishers, printers for printing logos on capsules, filling and packing equipment for gelatin capsules.

3) SPARE PARTS SUPPLY

We cooperate with manufacturers of automatic and semiautomatic capsule machines and we are engaged in the supply of spare parts for any model of capsule machines purchased in our company.



A handwritten signature in blue ink that reads "Roman Tsibulsky".

Owner and supervisor

ROMAN TSIBULSKY