

# MINIPRESS.RU

Supplier of pharmaceutical equipment in Russia  
Founded in 1999

RUSSIA, 115035, MOSCOW, 17 PIATNITSKAYA St.

+7(495)364-38-08

info@minipress.ru

www.minipress.ru

## Commercial proposal

**Compiled:** 15.02.2022  
**Valid:** 12 month  
**Author:** Roman Tsibulsky

### Automac capsule machine NJ-06

#### DESCRIPTION

Automac powder filling machine for hard gelatin capsules of size 000,00,0,1,2,3,4,5, A, B, C, D, E. Productivity 60-120 capsules per minute. Automac Powder Filling Machine for Solid Two-Component Gelatin Capsules. Optimum for use in the production of capsules. Weight 150 kg. Used for all standard sizes of hard gelatin capsules.

Made of stainless steel with GMP standard. We provide detailed instructions for setting up the capsule machine. Before sending the client capsule machine, machines are checked and tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. The price is indicated taking into account customs payments in Russia and delivery to the city of the client.

#### MANUFACTURER - CHINA



MINIPRESS.RU



MINIPRESS.RU

## SPECIFICATIONS

Automac capsule machine "NJ-06"

Material: stainless steel

Productivity: 60-120 capsules per minute

The sizes of gelatin capsules: 000,00,0,1,2,3,4,5, A, B, C, D, E.

Power supply: 0,5 kW, 220V, 50 Hz

Air pressure: 0,1 m<sup>3</sup> / minute 0,05 MPa

Overall dimensions: 605mm x 605mm x 1000mm

Weight: 150 kg

Weight with packing: 250 kg

A vacuum pump is included.

 [www.minipress.ru](http://www.minipress.ru)

## ADDITIONAL INFORMATION

Automac equipment for filling powder and granules of hard gelatin capsules. Empty capsules are poured into the hopper, then the machine guides the capsules in one position and fills the cells of the molds. Capsules are oriented downward with their lower parts. In the automac mode, the capsules are opened and their lower parts are filled with powder. The weight of the powder in all capsules is the same. After filling the capsule with powder, the upper and lower parts of the capsule are connected. Next, the ready-made capsules are poured into the tray. It is possible to download a video and view the entire cycle of the automac capsule filling machine. The operator loads the capsule halves into bins and powder, the machine works automatically. We recommend the use of gelatin capsules together with a polisher to remove powder particles from the surface of the capsules.

**PRODUCTION TIME**

**30-40 DAYS**

**PRICE**

**21186,44 USD**

**ESTIMATED DELIVERY TIME**

**30-40 DAYS**

**VAT**

**3813,56 USD**

**VALUE 25000.00 USD**

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank

### OUR SERVICES FOR CLIENTS:

#### 1) TRAINING OF COACHING TECHNOLOGY.

To all our customers who purchased equipment, we provide literature on the basics of encapsulation of powders. We carry out consultations by phone and by correspondence, we help to solve arising difficulties at manufacture on our equipment. We supply spare parts and consumables. We render services in the help on realization of the capsular equipment used.

#### 2) WE COMPLETE PRODUCTION CAPSULES.

We have a wide range of additional equipment for capsule production. Upon customer's request, we will select a model of capsule filling equipment, granulators, powder mixers, dust removal equipment and capsule polishers, printers for printing logos on capsules, filling and packing equipment for gelatin capsules.

#### 3) SPARE PARTS SUPPLY

We cooperate with manufacturers of automac and semiautomac capsulators and we are engaged in the supply of spare parts for any model of capsule machines purchased in our company.



*Roman Tsibulsky*

Owner and supervisor

**ROMAN TSIBULSKY**