

MINIPRESS.RU

Supplier of pharmaceutical equipment in Russia
Founded in 1999

RUSSIA, 115035, MOSCOW, 17 PIATNITSKAYA St.

+7(495)364-38-08

info@minipress.ru

www.minipress.ru

Commercial proposal

Compiled:	15.02.2022
Valid:	12 month
Author:	Roman Tsibulsky



Automatic capsule machine NJ-120

DESCRIPTION

Automatic machine for filling powder with hard gelatin capsules of size 00,0,1,2,3,4. Maximum productivity for capsules with granules 55000 capsules per hour. Maximum productivity for capsules with powder 72000 capsules per hour. Automatic Powder Filling Machine for Solid Two-Component Gelatin Capsules. Optimum for use in the production of capsules. Weight 1150 kg. Used for all standard sizes of hard gelatin capsules.

Maximum automation of the process. Made of stainless steel with GMP standard. We provide detailed instructions for setting up the capsule machine. Before sending the client capsule machine, machines are checked and tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. The price is indicated taking into account customs payments in Russia and delivery to the city of the client.

MANUFACTURER - CHINA



SPECIFICATIONS

Automac capsule machine "NJ-120"

Number of stations: 10

Cells in one station: 9

Material: stainless steel

Productivity for powders: 72000 capsules per hour

Productivity for granules: 55000 capsules per hour

The sizes of gelatinous capsules: 00,0,1,2,3,4

Power supply: 5.2 kW, 380 / 220V, 50 Hz

Air pressure: 0.1 m³ / minute, 0.05 MPa

Overall dimensions: 860mm x 940mm x 1900mm

Weight: 1150 kg

Weight with packing: 1300 kg

A vacuum pump is included.

ADDITIONAL INFORMATION

The principle of operation of automac capsule machines is as follows: hard gelatin capsules are collected in the hopper. According to the vacuum conveyor system, a group of capsules enters the corresponding cells of the turntable. There are processes of dividing the capsule into 2 parts, filling with powder, closing capsules. The setup process is rather complicated. Such mechanisms attach great importance to the quality of the capsules used. Repeatedly, we came across a situation where hard gelatin capsules, if stored incorrectly, or when a client bought a deliberately defective batch of capsules, led to a stoppage of the machine. This is due to the fact that the size and shape of the capsules may be impaired. The automac machine works at high speed and it is very difficult to track the penetration of substandard capsules. We recommend our customers purchasing a capsule automac machine to find the optimal quality gelatin capsules, and to observe all necessary measures for the storage of capsules. Maintenance of capsule automac machines requires the presence of personnel with appropriate skills. From the factory, the machines are delivered in a fully tuned state. When ordering a capsule machine, remember that the equipment for one type of capsule is included in the cost. Additional accessories for other sizes of capsules - for an additional payment. The location and the measuring plate are designed as a single unit to eliminate any deviations and avoid friction between the plate and the axis and thus improve filling accuracy, and also significantly extend the machine's performance. Simplicity and convenience of dismantling, installation and cleaning, various models of molds can replace each other, molds can be mutually replaced on the same machine to meet different performance requirements. The dust collector and vacuum tube, as well as the exhaust pipe, are installed inside the machine to avoid clogging the duct, breakage, leakage and other similar phenomena. The machine is more convenient for cleaning the platform. We recommend the use of gelatin capsules together with a polisher to remove powder particles from the surface of the capsules.

PRODUCTION TIME **30-40 DAYS**
ESTIMATED DELIVERY TIME **30-40 DAYS**

PRICE **47457,63 USD**
VAT **8542,37 USD**

VALUE 56000.00 USD

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank

OUR SERVICES FOR CLIENTS:

1) TRAINING OF COACHING TECHNOLOGY.

To all our customers who purchased equipment, we provide literature on the basics of encapsulation of powders. We carry out consultations by phone and by correspondence, we help to solve arising difficulties at manufacture on our equipment. We supply spare parts and consumables. We render services in the help on realization of the capsular equipment used.

2) WE COMPLETE PRODUCTION CAPSULES.

We have a wide range of additional equipment for capsule production. Upon customer's request, we will select a model of capsule filling equipment, granulators, powder mixers, dust removal equipment and capsule polishers, printers for printing logos on capsules, filling and packing equipment for gelatin capsules.

3) SPARE PARTS SUPPLY

We cooperate with manufacturers of automatic and semiautomatic capsule machines and we are engaged in the supply of spare parts for any model of capsule machines purchased in our company.



Owner and supervisor

ROMAN TSIBULSKY