

# MINIPRESS.RU

Supplier of pharmaceutical equipment in Russia  
Founded in 1999

RUSSIA, 115035, MOSCOW, 17 PIATNITSKAYA ST.

+7(495)364-38-08

info@minipress.ru

www.minipress.ru

## Commercial proposal

<b>Compiled:</b>	15.02.2022
<b>Valid :</b>	12 month
<b>Author:</b>	Roman Tsibulsky

### Mill for non-standard tablets RW-06

#### DESCRIPTION

Laboratory mill with continuous loading of 0-grade tablets. The design allows you to control the fraction of the resulting material at the output. Grinding with control of the powder fraction. Productivity from 20 to 50 kg per hour. Rotational speed of hammers is 1000 rpm. Portable, economical and easy to operate the machine is used in the medical and food industry for powder production, in the chemical industry and other industries. We provide a full range of services: installation, training, start-up, repair.

Simple mechanical design. Instruction for use in English and Russian. Complies with the GMP standard. We provide detailed instructions for setting up this model of shredders. Before sending the machine for grinding raw materials are tested and tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. Price is indicated with the delivery to the city of the client.

#### MANUFACTURER - CHINA



MINIPRESS.RU

MINIPRESS.RU



+7 495 505 10 10

## SPECIFICATIONS

Mill for non-standard tablets "RW-06"  
 Feeding material into the grinding zone: continuous  
 Hammer rotation speed: 1000 rpm  
 Air cooling: room temperature  
 The resulting powder fraction: 1000 to 200 microns  
 Changeable disc sieves: 3-5 pieces  
 Electric motor: 220V, 50Hz  
 Power: 1,5 kW  
 Dimensions: 470mm x 700mm x 800mm  
 Weight: 60 kg  
 Weight with packing: 80 kg

 [www.minipress.ru](http://www.minipress.ru)

## ADDITIONAL INFORMATION

This shredder is suitable for grinding o-grade pellets, to obtain a granular mass. A large selection of perforated conical nets with different cell sizes is included in the delivery. Customers can replace cones on their own to get the right size powdered powder fraction. This design has no competition for today, allowing to control the fraction, since in some cases, grinding into dust does not give the necessary effect for the further technological chain. In the kit there are also several fuses that protect against voltage fluctuations and a brush for cleaning the inner and outer surfaces of the chopper.

**PRODUCTION TIME** 30-40 DAYS  
**ESTIMATED DELIVERY TIME** 30-40 DAYS

**PRICE** 6779,66 USD  
**VAT** 1220,34 USD

**VALUE 8000.00 USD**

### OUR SERVICES FOR CLIENTS:

- 1) **TRAINING OF RAW MATERIAL TECHNOLOGIES.**  
 We can train you on our equipment before buying. Testing of materials.  
 We provide consultations on organization of raw materials grinding, with a full range of additional equipment and consumables.
- 2) **WE COMPLETE PHARMACEUTICAL AND CONFECTIONERY PRODUCTION.**  
 In Russia we have large and desktop mills that we provide to our customers. You can always contact us with your ideas and requests.
- 3) **DELIVERY OF EQUIPMENT FOR RAW MATERIAL.**  
 The products fully comply with the new GMP requirements. The mill output is accompanied by strict quality control, verification of all documentation, well-organized after-sales service of the equipment.
- 4) **SPARE PARTS SUPPLY**  
 The manufacturer can also offer a model of mills and grinders, based on the goals and objectives of the customer. We will assist in the maintenance and repair of any equipment for grinding powders. Service center in Russia for the CIS countries.

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank



*Roman Tsibulsky*

Owner and supervisor

**ROMAN TSIBULSKY**