

# MINIPRESS.RU

Supplier of pharmaceutical equipment in Russia  
Founded in 1999

RUSSIA, 115035, MOSCOW, 17 PIATNITSKAYA ST.

+7(495)364-38-08

info@minipress.ru

www.minipress.ru

## Commercial proposal

Compiled:	15.02.2022
Valid:	12 month
Author:	Roman Tsibulsky


### V-shaped mixer for powders VM-50

#### DESCRIPTION

Highly effective mixer for mixing dry powder material in the pharmaceutical, food, chemical and other industries. Thanks to the unique design, uniform mixing occurs. Useful volume - 50 liters. Mixing capacity is made of stainless steel with polished inner and outer walls. Implemented by GMP.

We provide detailed instructions for setting up a mixer of powders. Before shipment to the customer, the mixer is tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. Price is indicated with the delivery to the city of the client.

#### MANUFACTURER - CHINA



MINIPRESS.RU



MINIPRESS.RU



+7 495 505 10 10

## SPECIFICATIONS

Mixer V-shaped for powders "VM-50"  
 Type: non-contact  
 Application: mixing of dry powders  
 Volume: 50 liters  
 Possible working load: 20 liters  
 Optimum in quality of loading: 14 liters  
 Rotation speed: 15 rpm  
 Coefficient of loading by powder: 40-60%  
 Mixing time: 15-20 minutes  
 Power: 0,75 kW  
 Power supply: 220V or 380V at customer's choice  
 Overall dimensions: 1640mm x 670mm x 1720mm  
 Weight: 180 kg  
 Shipping Weight: 270 kg

 [www.minipress.ru](http://www.minipress.ru)

## ADDITIONAL INFORMATION

The mixer is made of stainless steel. The V-shaped high-efficiency mixer has no blind areas in the cylinders at high speed, and with minimal mixing time. The body, internal and external walls of polished stainless steel, have an elegant appearance. Loading and unloading of the powder can be done either manually or by means of a vacuum conveyor. To unload the powder, there are quick-release valves. This mixer is used for mixing powdered materials in the pharmaceutical and other industries. From the electric motor, the shaft should be set in motion. The best industrial device for mixing dry free-flowing powders.

**PRODUCTION TIME** 30-40 DAYS  
**ESTIMATED DELIVERY TIME** 30-40 DAYS

**PRICE** 7203,39 USD  
**VAT** 1296,61 USD

**VALUE 8500.00 USD**

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank

### OUR SERVICES FOR CLIENTS:

#### 1) TRAINING OF MIXING TECHNOLOGY.

To all our customers who purchased equipment, we provide literature on the basics of mixing. We carry out consultations by phone and by correspondence, we help to solve arising difficulties at manufacture on our equipment. We supply spare parts and consumables. We render services in the help on realization of the capsular equipment used.

#### 2) WE COMPLETE PRODUCTION CAPSULES.

We have a wide range of additional equipment for mixing production. Upon customer's request, we will select a model of capsule filling equipment, granulators, powder mixers, dust removal equipment and capsule polishers, printers for printing logos on capsules, filling and packing equipment for gelatin capsules.

#### 3) SPARE PARTS SUPPLY

We cooperate with manufacturers of mixers and we are engaged in the supply of spare parts for any model of mixers purchased in our company.



*Roman Tsibulsky*

Owner and supervisor

**ROMAN TSIBULSKY**