

MINIPRESS.RU

Supplier of pharmaceutical equipment in Russia
Founded in 1999

RUSSIA, 115035, MOSCOW, 17 PIATNITSKAYA St.

+7(495)364-38-08
info@minipress.ru
www.minipress.ru

Commercial proposal

Compiled:	15.02.2022
Valid:	12 month
Author:	Roman Tsubulsky


Automac Filling and Sealing Machine for Sterile Ampoules ABF-4D-1

DESCRIPTION

Automac machine for filling liquids and sealing glass sterile ampoules of a closed type for pharmaceutical production. The productivity is 3600-8100 ampoules per hour. Dosage of ampoules in several versions for 1-2 ml, 5-10 ml, 20 ml, 1-10 ml, 5-20 ml. Weight 500 kg. This model uses peristaltic pumps and an electronic control system with a LCD screen. Automac, economical and easy to operate the machine is used in the medical, food and chemical industries. We provide a full range of services: installation, training, start-up, repair.

Instruction for use in English and Russian. Complies with the GMP standard. We provide detailed instructions for setting up this model of the machine for filling and sealing ampoules. Before sending the machine for sealing, machines are tested and tested in production. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. The price is indicated taking into account customs payments in Russia and delivery to the city of the client.

MANUFACTURER - CHINA



SPECIFICATIONS

Automac filling and sealing machine for sterile ampoules "ABF-4D-1"

Modifications: 1-2 ml, 5-10 ml, 20 ml, 1-10 ml, 5-20 ml

Productivity (1-2 ml): 7200 - 8100 ampoules per hour.

Productivity (5-10 ml): 5400 - 6000 ampoules per hour.

Productivity (20 ml): 3600 ampoules per hour.

Productivity (1-10 ml): 5400 - 7200 ampoules per hour.

Productivity (5-20 ml): 5400 - 600 ampoules per hour.

Power: 0.55 kW

Gas consumption: 2-2.5 m³ per hour

Air consumption: 0,08-0,1 MPa

Supporting combustion gas consumption: 1 kg per hour

Oxygen consumption: 7 kg per hour

Power supply: 220V or 110V,

Dimensions: 1980mm x 1000mm x 1500mm

Weight: 500 kg

Weight with packing: 600 kg

 www.minipress.ru

ADDITIONAL INFORMATION

Automac machines for boling and sealing sterile ampoules of closed type. Models with the ABF-4D-1 marking are supplied with glass pumps with metering devices and a push-buon control system. There are more expensive ABF-4D-2 models in which one peristalc pump is installed for each of the four ampoules and an electronic control system with a LCD display. Otherwise, all models dier from each other with minor opons, dosage and performance. The fence of the closed ampoules comes from the hopper, then they enter the unloading zone through the transporter. The ame melts the p of the ampoule and boles are produced. The operator adjusts the dosage, within the capabilities of each modicaon of the machine. Filling in ampoules is carried out by a group of 4 ampoules. The machine brings a group of ampoules under the needles, the conveyor stops and the needles fall into ampoules. Next comes the lling of ampoules. Then the glass ampoules are moved to the burner zone. The combuson temperature and the intensity of the ame are also regulated. Aer the neck of the ampoule has been heated by special clamps, the apex of the ampoule is formed. The p has a rounded contour of regular uniform shape. The height of all ready-made ampoules and roundness is repeated constantly, without disncon. Attenon, the machine for lling and sealing sterile glass ampoules does not require a machine for washing and drying ampoules. It is necessary to use closed glass ampoules.

PRODUCTION TIME

30-40 DAYS

ESTIMATED DELIVERY TIME

30-40 DAYS

PRICE

18644,07 USD

VAT

3355,93 USD

VALUE 22000.00 USD

Payment is made in rubles to the account in Moscow at the rate of the RF Cental Bank

OUR SERVICES FOR CLIENTS:

1) TRAINING TECHNOLOGIES OF BOTTLING AND BATCHING AMPUL.

We can train you on our equipment before buying. Testing of materials. We provide consultations on filling and sealing of ampoules, with a full range of additional equipment and consumables.

2) WE COMPLETE PHARMACEUTICAL AND CONFECTIONERY PRODUCTION.

We have automatic and desktop machines for filling and sealing ampoules in Russia, which we provide to our customers. You can always contact us with your ideas and requests.

3) SUPPLY OF EQUIPMENT FOR BOTTLING AND STORING OF GLASS AMPOULES.

The products fully comply with the new GMP requirements. The release of equipment for sealing glass ampoules is accompanied by strict quality control, verification of all documentation, well-organized after-sales service of equipment.

4) SPARE PARTS SUPPLY

The manufacturer can also offer models of automatic machines for filling and sealing ampoules, based on the goals and objectives of the customer. We will assist in the maintenance and repair of any equipment for sealing ampoules. Service center in Russia for the CIS countries.



Roman Tsibulsky

Owner and supervisor

ROMAN TSIBULSKY