

# Commercial proposal

Compiled:

15.02.2022

Valid:

12 month

Author:

Roman Tsibulsky

## Table-top blister machine ITA-03

### DESCRIPTION

Table blister packing machine, to fill the finished blisters with tablets and capsules, followed by sealing with foil. The equipment has the function of embossing the expiry date and the production date. Suitable for laboratory use, as well as for small production of blisters. It is made of stainless steel. Productivity 500-600 blisters per hour.

Producon Italy. We provide detailed instructions for setting up a blister machine. Before sending the customer a blister machine is tested in producon. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. Price is indicated with the delivery to the city of the client.

### MANUFACTURER-ITALY

**MINIPRESS.RU****MINIPRESS.RU**



## SPECIFICATIONS

Machine for welding foil to PVC "ITA-03"

The equipment allows you to seal the blisters of standard shapes.

Used forms must be pre-ordered by sending a design.

The machine is completely made of stainless steel and aluminum.

Equipped with an electronic thermoregulator, which automatically regulates the temperature of the foil.

Allows you to manually adjust the degree of heat of phase welding, depending on the thickness of the material used.

Productivity: 500-600 blisters per hour

Overall dimensions: 300 x 260 x 460 mm

Power supply: 220 V, 300 W, 50/60 Hz.

Shipping Weight: 10 kg

Blister preparations can be purchased from us.

[www.minipress.ru](http://www.minipress.ru)

## ADDITIONAL INFORMATION

The machine is semi-automatic, it is installed on the table. Before use, the customer must order from the manufacturer an aluminum matrix according to his design, and have ready blanks of empty blisters. Aluminum foil is partially welded to the blisters. Manually it is necessary to place an empty blister in an aluminum mold, fill all the blister cells with tablets or capsules, cover with aluminum foil from above. Move the aluminum mold inside, turn on the temperature regulator and set the desired heating temperature, then manually press the foil against the PVC. After a pause, you need to remove the aluminum mold and get the finished blister.

PRODUCTION TIME	30-40 DAYS	PRICE	3389,83 USD
ESTIMATED DELIVERY TIME	30-40 DAYS	VAT	610,17 USD

**VALUE 4000.00 USD**

### OUR SERVICES FOR CLIENTS:

Payment is made in rubles to the account in Moscow at the rate of the RF Central Bank

#### 1) TRAINING OF BLISTERING TECHNOLOGY.

To all our customers who purchased equipment, we provide literature on the basics of packing tablets and capsules in blisters. We carry out consultations by phone and by correspondence, we help to solve arising difficulties at manufacture on our equipment. We supply spare parts and consumables. We render services in the help on realization of the blister equipment of the user.

#### 2) WE COMPLETE MANUFACTURE OF BLISTERS.

We have a wide range of additional equipment for blister production. At the request of the client we will select the model of an automatic blister machine, daters and numbering machines, packing and packing equipment for gelatin capsules and tablets.

#### 3) SPARE PARTS SUPPLY

We cooperate with manufacturers of automatic and semi-automatic blister machines and we are engaged in the supply of spare parts for any model of blister equipment purchased in our company.



*Roman Tsibulsky*

Owner and supervisor

**ROMAN TSIBULSKY**