

Commercial proposal

Compiled:	15.02.2022
Valid:	12 month
Author:	Roman Tsibulsky


Thickness analyzer of the walls of gelatin capsules HD-01

DESCRIPTION

The thickness analyzer is used to measure the thickness of tablets or gelatin capsules. The analyzer measures the thickness of the top and the wall of the capsule. The device gives the most accurate result. High-grade cengrade scale. We provide detailed instrucons for seng up this model of the thickness analyzer.

Before sending, the analyzers are tested in producon. Completeness and working capacity is guaranteed. We maintain a stock of parts and consumables in the warehouse. The price is indicated taking into account customs payments in Russia and delivery to the city of the client.

MANUFACTURER - CHINA



MINIPRESS.RU

MINIPRESS.RU

SPECIFICATIONS

Analyzer thickness of the walls of gelan capsules "HD-01"

Test range: (0-5) mm

Accuracy of the test: ± 0.01 mm

Type of capsule: 00,0,1,2,3,4,5

Size with packing: 260mm x 140mm x 140mm

Shipping Weight: 1.5 kg

 www.minipress.ru

ADDITIONAL INFORMATION

The analyzer measures the thickness of the top and the walls of a hard gelan capsule. The hand-held device gives the most accurate measurement result. High-grade cengrade scale. Suitable for hard gelan capsules of sizes 000,00,0,1,2,3,4,5.

PRODUCTION TIME 30-40 DAYS

ESTIMATED DELIVERY TIME 30-40 DAYS

PRICE 1101,69 USD

VAT 198,31 USD

VALUE 1300.00 USD

Payment is made in rubles to the account in Moscow at the rate of the RF Cental Bank

OUR SERVICES FOR CLIENTS:

1) TRAINING OF MEDICINE VERIFICATION TECHNOLOGY.

To all our customers who purchased equipment, we provide literature on the basics of encapsulaon of powders. We carry out consultaons by phone and by correspondence, we help to solve arising dicules at manufacture on our equipment. We supply spare parts and consumables. We render services in the help on realizaon of the capsular equipment used.

2) WE COMPLETE PRODUCTION CAPSULES.

We have a wide range of additional equipment for capsule producon. Upon customer's request, we will select a model of capsule lling equipment, granulators, powder mixers, dust removal equipment and capsule polishers, printers for prinng logos on capsules, lling and packing equipment for gelan capsules.

3) SPARE PARTS SUPPLY

We cooperate with manufacturers of automac and semiautomac capsulators and we are engaged in the supply of spare parts for any model of capsule machines purchased in our company.



Roman Tsibulsky

Owner and supervisor

ROMAN TSIBULSKY